




AlphaLab SnR/Neuro Omega

AO_MatlabTool

User's Manual

Version 1.2



ALPHA OMEGA

Home Office: P.O. Box 810, Nazareth Illit 17105, Israel

Email: info@alphaomega-eng.com

Web site: www.alphaomega-eng.com

Copyright □ 2010-2013 Alpha Omega Engineering LTD

Proprietary Notice

All Rights Reserved

This publication, or parts thereof, contains information proprietary to Alpha Omega Engineering and may not be reproduced, recorded or transmitted in any form, by any method, for any purpose without written permission from Alpha Omega Engineering.

Alpha Omega Engineering LTD. makes no warranty, expressed or implied, including but not limited to any implied warranties of merchantability or fitness for a particular purpose, regarding these materials and makes such materials available solely on an “as-is” basis.

In no event shall Alpha Omega Engineering LTD. be liable to anyone for specific, collateral, incidental, or consequential damages in connections with or arising from purchase or use of these materials. The sole and exclusive liability to Alpha Omega Engineering, regardless of the form of actions, shall not exceed the purchase price of the materials described herein.

AlphaLab SnR, Neuro Omega, and Alpha Omega are trademarks of Alpha Omega Engineering. Any other product and/or company names are trademarks of their respective owners. Due to continuing product development, specifications, capabilities and features are subject to change without notice.

Contact Information:

ALPHA OMEGA ENGINEERING

P.O. Box 810

Nazareth Illit 17105

Israel

Tel: +972-4-656-3327

Fax: +972-4-657-4075

info@alphaomega-eng.com

support@alphaomega-eng.com

<http://www.alphaomega-eng.com>

U.S. Office: Toll Free: 1-877-919-6288, Fax: 1-877-471-2055

Europe Office: Toll Free: 00-800-2-574-2111, Tel +49-7-251-440-6620

Table of contents

| | |
|--|----|
| 1. Introduction | 1 |
| 1.1. About This Manual..... | 1 |
| 1.2. Reference Documents..... | 1 |
| 1.3. List of Abbreviations | 1 |
| 1.4. Matlab Tool Overview | 1 |
| 1.5. Matlab Tool Command list..... | 1 |
| 2. Getting started..... | 3 |
| 3. Matlab Functions and Usage | 5 |
| 4. Appendix-1 - Matlab function return cases | 19 |
| 5. Appendix-2: AlphaLab SnR\Neuro Omega channel list | 20 |

1. Introduction

1.1. About This Manual

This manual describes the features and usage of the Matlab data access and control tool for AlphaLab SnR\Neuro Omega system.

This manual will use the AlphaLab SnR when describing the function, but all described functions also apply to Neuro Omega System.

1.2. Reference Documents

AlphaLab SnR Manual
Neuro Omega Manual

1.3. List of Abbreviations

SnR- Stimulation and Recording

1.4. Matlab Tool Overview

Matlab Tool is a set of commands allowing multi user communication with the AlphaLab SnR\Neuro Omega system. Each command is a Matlab function that allows the user to control certain system features.

Matlab Tool is compatible with Windows XP/7 32 bit.

1.5. Matlab Tool Command list

The table below shows a list of available Matlab function and a brief description of each one. The detailed usage of each command is explained later in the manual.

| <i>Function</i> | <i>Description</i> |
|------------------------|--|
| AO_StartConnection | Connect Matlab to AlphaLab SnR |
| AO_IsConnected | Check if Matlab is connected to the AlphaLab SnR |
| AO_CloseConnection | Close connection between Matlab and the AlphaLab SnR |
| AO_AddBufferingChannel | Acquire data for the specified channel |
| AO_GetAlignedData | Get aligned data from several channels |

| <i>Function</i> | <i>Description</i> |
|-----------------------------|---|
| AO_GetChannelData | Get data for the specified channel |
| AO_ClearChannelData | Clear all buffered data |
| AO_GetNextBlock | Get next block of data |
| AO_SendBlock | Send stream format data |
| AO_StartSave | Start saving mpx file on the AlphaLab SnR |
| AO_SetSaveFileName | Set the mpx file name saved by AlphaLab SnR |
| AO_SetSavePath | Set the path of the directory where the lsx/mpx files are saved by AlphaLab SnR |
| AO_StopSave | Stop saving on the AlphaLab SnR |
| AO_SetStimulationParameters | Set stimulation parameters |
| AO_StartStimulation | Start stimulation |
| AO_StopStimulation | Stop stimulation |
| AO_SendDout | Send digital output to a specific Port in AlphaLab SnR |

2. Getting started

This section describes the preparation work that needs to be done in order to be able to use the Matlab tool. Follow the steps below:

1. Install Matlab Tool
 - a. Run the supplied setup file, and follow instructions
2. Start Matlab
 - a. Open Matlab (Under Windows 7, you may need to run Matlab as Administrator, or change the user settings to lower Parental control)
 - b. Set the working directory path in Matlab to the installed Matlab Tool Directory (Usually C:\Program Files\AlphaOmega\AO_MatlabTool)
3. Set up compiler and compile the mex file as follow:
 - a. In Matlab command window type “**mex -setup**” and press enter; A Matlab message will appear in the command window:

‘Would you like mex to locate installed compilers [y]/n?’
 - b. Press **n**
The Matlab will suggest a list of all compilers it supports
 - c. Choose a version of Microsoft visual such as Microsoft free visual c++ 2008 or 2010. *Note: If you do not have it on your PC, you need to install it before continuing (express mode is downloaded for free)*
 - d. Continue procedure of choosing the compiler by answering the questions. This is for the path validation, so answer “y” to all questions if the path is correct.
4. Compile the Mex files as follows:
 - a. In Matlab command window type “**mex MexFileEthernetStandAlone.cpp Include\ethernetStandAlone.lib**” and press enter; this command will compile the MexFileEthernetStandAlone.cpp and create a mex file called MexFileEthernetStandAlone.mexw32, and this completes the installation.

5. Test the installation as follow:

- a. In Matlab command window type: “**AO_IsConnected**“ and press enter. If no compilation error, which is usually indicated by red colored messages, you are all set and ready to go.

Otherwise,

- b. In Matlab command window type
! dependency walker\ depends.exe MexFileEthernetStandAlone.mexw32
and press enter, and this should indicate which dll files are missing in order to run MexFileEthernetStandAlone.mexw32. If you are unable to obtain the missing dll's, contact Alpha Omega support.

Note:


In case of an old Matlab version, a compiling error will appear while compiling the cpp file. In this case:

- a. Using any text editor, edit ‘MexFileEthernetStandAlone.cpp’ file
- b. Go to the line `#define NEW_MATLAB`
- c. Add `//` at the beginning of the line to comment it out. (`///#define NEW_MATLAB)`

3. Matlab Functions and Usage

The table below gives complete syntax, description, and an example of each Matlab function.

| Function | Function syntax and example |
|---------------------------|--|
| AO_StartConnection | <p>[Result] = AO_StartConnection(DspMac,PcMac,AdpaterIndex)</p> <p>Used to connect Matlab to AlphaLab SnR system</p> <p><u>Function parameters:</u></p> <p>DspMac: String of 6 hex values. This is the mac address of the SnR system</p> <p>PcMac: String of 6 hex values. This is the mac address of matlab PC. This value must be different from the DspMac and the Mac address of the PC running the SnR software, it is recommended to find the actual mac address of the Matlab computer and use it.</p> <p>AdpaterIndex: Constant int value of -1</p> <p>Result: Function return is an integer, 0 = no function errors, other number indicate function error (See appendix 1)</p> <p>It is preferable to ensure connection was done successfully by calling the function AO_IsConnected</p> <hr/> <p>%Example</p> <pre>PCMac='00:01:02:2:BD:BD'; DSPMac='00:21:ba:07:ab:9e'; Activeadapter=-1; retStartConnection=AO_startConnection(DSPMac,PCMac,Activeadap ter);</pre> <p>The following code should be added to insure proper connection:</p> |

| Function | Function syntax and example |
|-----------------------|--|
| | <pre> for j=1:100, pause(1); ret=AO_IsConnected; if ret==1 'The System is Connected' break; end end </pre> <hr/> <p>If connection was done successfully , the PcMac address will appear in AlphaLab SnR GUI => Help => User info</p>  |
| AO_IsConnected | <p>[Results] = AO_IsConnected()</p> <p>This function checks if Matlab is connected to the AlphaLab SnR</p> <p>Returns 1 if the system is connected, otherwise returns 0</p> <hr/> <p>%Example</p> <pre> for j=1:100, </pre> |

| Function | Function syntax and example |
|-------------------------------|---|
| | <pre> pause(1); ret=AO_IsConnected; if ret==1 'The System is Connected' break; end end </pre> |
| AO_CloseConnection | <p>[Result] = AO_CloseConnection()</p> <p>Function used to close connection between Matlab and the AlphaLab SnR system</p> <p>Result: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix 1)</p> <hr/> <p>% Example</p> <pre> Result=AO_CloseConnection(); if (Result==0) display('Connection closed successfully'); else display('Connection close error'); end </pre> |
| AO_AddBufferingChannel | <p>[Result] = AO_AddBufferingChannel(ChannelID,SamplingRate)</p> <p>Function used to gather data for the channel defined in ChannelID</p> <p><u>Function parameters:</u></p> <p>ChannelID: The channel ID we want to gather data for</p> <p>SamplingRate: Sampling rate of the channel in samples/sec (this used strictly to calculate buffering information).</p> <p>Result: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix 1)</p> |

| Function | Function syntax and example |
|--------------------------|---|
| | <p>Notes:</p> <ul style="list-style-type: none"> • The size of buffering is 15sec* SamplingRate • The function stores the data using First In First Out (FIFO) mechanism • The data value is A/D value including gain <hr/> <p>%Example</p> <pre>ChannelID=10256; % set the channel number SamplingRate=44642.8571428; %set the sampling rate of the channel AO_AddBufferingChannel(ChannelID,SamplingRate)% start gathering data for channel 10256 with sampling rate SamplingRate</pre> |
| AO_GetAlignedData | <p>[Result,pData>DataCapture,TS_FirstSample]=AO_GetAlignedData(ChannelIdArr,ChannelCount)</p> <p>Function used to get aligned data for several channels.</p> <p><u>Function parameters:</u></p> <p>ChannelIdArr: Array of channels which we need to get data for, all channels must have the same sampling rate</p> <p>ChannelCount: Number of channels listed in ChannelIdArr</p> <p>Result: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix 1)</p> <p>pData: Array of data samples of all listed channels, be arranged in single row</p> <p>DataCapture: The amount of the useable data captured in the pData array</p> <p>TS_FirstSample:The timestamp of the first sample for each channel</p> <p>Notes:</p> <ul style="list-style-type: none"> • In order to get data you need to use AO_AddBufferingChannel first. • In order to get real time data u need to clear the buffered data first using AO_ClearChannelData function before using both commands |

| Function | Function syntax and example |
|--------------------------|--|
| | <p>AO_GetAlignedData, AO_GetChannelData otherwise u will get stored data</p> <ul style="list-style-type: none"> pData will contain samples of data for all channels ,the number of valid samples in this array is DataCapture so make sure that you only get DataCapture samples. Hence, number of samples for each channel is DataCapture divided by ChannelCount. In pData samples are arranged in a single array for all channels, starting with samples of the first channel listed in the ChannelIdArr , followed by other channels consecutively and in the same order. The data value is A\D including gain <hr/> <p>%Example</p> <pre>ChannelIdArr=[10000,10001,10002]; ChannelCount=3; [Result,pData,DataCapture,TS_FirstSample] = AO_GetAlignedData(ChannelIdArr,ChannelCount);% get aligned data from channels:10000,10001,10002 save them in the array pData, the alignment is done by time stamp TS_FirstSample</pre> |
| AO_GetChannelData | <p>[Result,pData,DataCapture] = AO_GetChannelData(ChannelId)</p> <p>Function used to get data for the specified channel.</p> <p><u>Function parameters:</u></p> <p>ChannelId: The channel id which we want to get data for pData: Array of data DataCapture: The amount of the useable data in the array Result: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix 1)</p> <p>Notes:</p> <ul style="list-style-type: none"> pData will contain a block of data |


| Function | Function syntax and example |
|----------------------------|---|
| | <ul style="list-style-type: none"> The format of the data block is as follow: byte 1-2 : SizeOftheBlock in words (1 word =2Byte) including the samples in this block byte 3 BlockType (in our case alwayes will be 'd' or 100) byte 4 Not used byte 5-6 ChannelNumber the id of the channel this block belongs to byte 7 Unit number ,this value valid only for segmented data byte 8 Not used byte 9-12 TimeStamp of the first sample of the block you will have to reorder them [byte10 byte9 byte12 byte11] byte 13-14 OverFlowCount the over flow of the time stamp – Future use byte 15-16 First sample byte 17-18 Second sample etc.... In order to calculate the number of samples in this channel do the following HeaderSize=14bytes HeadrSizeWord=14bytes/2 samplescount=SizeOftheBlock-HeaderSizeWord = (SizeOftheBlock-14)/2 <hr/> <p>%Example</p> <p>[Result,pData,DataCapture]=AO_GetChannelData(10256);</p> |
| AO_ClearChannelData | <p>[Result] = AO_ClearChannelData()</p> <p>Function used to clear buffered data by command AO_AddBufferingChannel</p> <p>Result: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix 1)</p> |

| Function | Function syntax and example |
|-------------------------------|---|
| | <p>Note:</p> <ul style="list-style-type: none"> In order to get real time data you need to clear the buffered data first using AO_ClearChannelData function before using both commands AO_GetAlignedData , AO_GetChannelData otherwise u will get stored data <hr/> <p>%Example</p> <p>AO_ClearChannelData()</p> |
| <p>AO_GetNextBlock</p> | <p>[Result,arraydata,realDataSizeWords]=AO_GetNextBlock(sizeOfArrayWords)</p> <p>Function used to get the next new block data, the data should be parsed using StreamFormat.h file</p> <p><u>Function parameters:</u></p> <p>sizeOfArrayWords: The max size of data the array can contain</p> <p>arraydata: Pointer to an array to hold the new data ,the data contain stream format in order to parse the data you need some Knowledge in our stream Format</p> <p>realDataSizeWords: The count of the data copied to the array data</p> <p>Result: Function returns an integer, 0 = no function errors, other number indicate function error (appendix 1)</p> <p>Note:</p> <ul style="list-style-type: none"> StreamFormat.h file is saved in the include directory <hr/> <p>%Example</p> <p>realDataSizeWords=zeros(1,1); [res,arraydata,realDataSizeWords]=AO_GetNextBlock(50000);</p> |

| Function | Function syntax and example |
|---------------------|--|
| AO_SendBlock | <p>[Result] = AO_SendBlock(ArrayData)</p> <p>Function used to send stream format data to the AlphaLab SnR system</p> <p><u>Function parameters:</u></p> <p>ArrayDat: Contain the data which will be sent to AlphaLab SnR system</p> <p>Result: Function return is an integer, 0 = no function errors, other number indicate function error (See appendix 1)</p> <p>Notes:</p> <ul style="list-style-type: none"> • This function for advanced users only. • Stream format is explained in StreamFormat.h file <hr/> <p>%Example</p> <pre>ArrayData=[7,1,2,3,4,5,6]; AO_SendBlock(ArrayData)</pre> |
| AO_StartSave | <p>[Result] = AO_StartSave();</p> <p>Function used to start saving mpx file by the AlphaLab SnR system</p> <p><u>Function parameters:</u></p> <p>Result: Function returns an integer, 0 = no function errors, other number indicate function error (appendix 1)</p> <p>Notes:</p> <ul style="list-style-type: none"> • The mpx file saved will contain the channels listed on the Data logging Options in the SnR • The filename and the saving path could be set before saving using Matlab commands: AO_SetSaveFileName, and AO_SetSavePath. |

| Function | Function syntax and example |
|---------------------------|---|
| | <ul style="list-style-type: none"> • Or by the parameters defined in the data logging(default) • When saving start the saving button in the Alpha Lab SnR GUI turns to red. See AlphaLab SnR Manual <hr/> <p>%Example:</p> <p>[Result] = AO_StartSave()% start saving on the SnR/Neuro Omeg</p> |
| AO_SetSaveFileName | <p>[Result] = AO_SetSaveFileName(FileName)</p> <p>Function used to set mpx file name saved by AlphaLab SnR system</p> <p><u>Function parameters:</u></p> <p>FileName: Contains the file name. File name must be less than 30 chars</p> <p>Result: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix1)</p> <hr/> <p>%Example:</p> <p>fileName='TestFile'; AO_SetSaveFileName(fileName)%set the file name as TestFile</p> <p>AO_StartSave()% start saving, the file name will be testFile</p> |
| AO_SetSavePath | <p>[Result] = AO_SetSavePath(Path)</p> <p>Function used to set the path of the directory to save in the mpx file</p> <p><u>Function parameters:</u></p> <p>Path: Contain the path of the directory for saving the files</p> <p>Result: Function return is an integer, 0 = no function errors, other number indicate function error</p> |

| Function | Function syntax and example |
|------------------------------------|---|
| | <p>(See appendix 1)</p> <hr/> <p>%Example</p> <pre>path='c:\logging_data\';%the path of the directory to save in AO_SetSavePath(path)%set the path of the saving to 'c:\logging_data\' AO_StartSave()%start saving, the file will be saved at 'c:\logging_data\'</pre> |
| AO_StopSave | <p>[Result] = AO_StopSave()</p> <p>Function used to stop saving by AlphaLab Snr/Neuro Omega system</p> <p>Result: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix 1)</p> <hr/> <p>%Example</p> <pre>fileName='TestFile'; path='c:\logging_data\';;% AO_SetSaveFileName(fileName); AO_SetSavePath(path); AO_StartSave();% pause(100); AO_StopSave();</pre> <p>➔saving will be done for 100 sec, the mpx file name is TestFile in the 'c:\logging_data\'</p> |
| AO_SetStimulationParameters | <p>[Result] = AO_SetStimulationParameters (FirstPhaseDelay_mS,FirstPhaseAmpl_mA,FirstPhaseWidth_mS, SecondPhaseDelay_mS,SecondPhaseAmpl_mA,SecondPhaseWidth_mS,Freq_hZ,Duration_sec,ReturnChannel,channelnumber);</p> <p>Function used to set stimulation parameters</p> |

| Function | Function syntax and example |
|---|--|
|  | <p><u>Function parameters:</u></p> <p>FirstPhaseDelay_mS: First phase delay in mSec (1) FirstPhaseAmpl_mA: First phase amplitude (4) FirstPhaseWidth_mS : The width of the first phase (3) SecondPhaseDelay_mS: Second phase delay in mSec (2) SecondPhaseAmpl_mA: Second phase amplitude (6) SecondPhaseWidth_mS: Second phase width (5) Freq_hZ: The stimulation frequency Duration_sec: Duration of the stimulation after which stimulation stops ReturnChannel: The ID of the channel we want to return the stimulation with(set -1 for Global return) channelnumber: The channel we want to start stimulation on Results: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix1)</p> <p>Note:</p> <ul style="list-style-type: none"> This function should be called before starting stimulation, otherwise stimulation will be done using the parameters defined in the SW GUI <hr/> <p>%Example</p> <pre> FirstPhaseDelay_mS=1.1;%the delay of the first phase FirstPhaseAmpl_mA=-3.5;%the amp of the first phase FirstPhaseWidth_mS=0.5;%the width of the first phase SecondPhaseDelay_mS=1.5;%the delay of the second phase SecondPhaseAmpl_mA=1.5;%the amp of the second phase SecondPhaseWidth_mS=0.2;%the width of the second phase Freq_hZ=10;%the frequency of the stimulation Duration_sec=30;%duration of the stimulation ReturnChannel=10001;%the ID of the channel we want to return the stimulation with channelnumber=10000;%the channel we want to start stimulation in AO_SetStimulationParameters(FirstPhaseDelay_mS,FirstPhaseAmpl_ </pre> |

| Function | Function syntax and example |
|----------------------------|---|
| | mA,FirstPhaseWidth_mS,SecondPhaseDelay_mS,SecondPhaseAmpl_mA,SecondPhaseWidth_mS,Freq_hZ,Duration_sec,ReturnChannel,channelnumber);%set stimulation params |
| AO_StartStimulation | <p>[Result] = AO_StartStimulation(ChannelNumber);</p> <p>Function used to start stimulation using AlphaLab SnR system</p> <p><u>Function parameters:</u></p> <p>ChannelNumber: The ID of the channel used for stimulation</p> <p>Result: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix 1)</p> <p>Notes:</p> <ul style="list-style-type: none"> • The stimulation parameters should be set before stimulation using Matlab command AO_SetStimulationParameters, otherwise it will use the SW UI stimulation parameters • Be aware when stimulation is done with more than one channel, that the set stimulation parameters refer to the stimulation source • When stimulation is on, stimulation button on the GUI turn to red <hr/> <p>%Example</p> <p>ChannelNumber=10000; AO_StartStimulation(ChannelNumber);</p> |
| AO_StopStimulation | <p>[Results]=AO_StopStimulation(ChannelNumber);</p> <p>Function used to stops stimulation to the stimulation source of the ChannelNumber</p> <p><u>Function parameters:</u></p> |

| Function | Function syntax and example |
|-------------|---|
| | <p>ChannelNumber: The ID of the channel used for stimulation</p> <p>Results: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix1)</p> <p>Note: In order to stop stimulation in all channels use: ChannelNumber= -1</p> <hr/> <p>%Example</p> <pre>ChannelNumber=10000; AO_StartStimulation(ChannelNumber);</pre> |
| AO_SendDout | <p>[Results] =</p> <p>AO_SendDout(DigitalChannelNumber,mask,Value)</p> <p>Function used to send Dout to a specific Port</p> <p><u>Function parameters:</u></p> <p>DigitalChannelNumber: The ID of the digital port</p> <p>Mask: Is an 8 bit hex number input as a string. It indicates which bits are to be changed. Any bit set to 1 will be changed to the corresponding bit in the value, and any bit set to 0 will not be changed.</p> <p>Value: Can be any number between 0 and 2^8-1.</p> <p>Results: Function returns an integer, 0 = no function errors, other number indicate function error (See appendix1)</p> <hr/> <p>%Example</p> <pre>DigitalChannelNumber=11701 ; %channel ID mask='0xFF'; %the mask (11111111) value=0; %the value (00000000)</pre> <p>Result = AO_SendDOut(DigitalChannelNumber,mask,value); %Initialize all bits of port 11701 to 0</p> |

| Function | Function syntax and example |
|----------|--|
| | <pre> mask='0x05'; %the mask (00000101) value=3; %the value (00000011) Result = AO_SendDOut(DigitalChannelNumber,mask,value); %set port 11701 %====> The output of the bits on port 11701 will be '0000 0001' Mask=00000101 Value=00000011 Port= 00000001 </pre> |

4. Appendix-1 - Matlab function return cases

| Function return | Result |
|-----------------|---|
| 0 | No compiling error |
| 1 | The frequency must be bigger than 20 |
| 2 | Duration of the signal is longer than 70 msec and its must be shorter |
| 3 | Null parameter |
| 4 | dll is not connected |

5. Appendix-2: AlphaLab SnR\Neuro Omega channel list

Base addresses for each channel is defined in the workspace for example it could start with 10000 for first system, 20000 for second system, 30000 for third system and so on. Therefore, when addressing a specific channel you should address 10000+ channel number.

| Number | | | | |
|--|-----|---|-----------------------------------|--------------|
| From | To | Signal | From / To | Type |
| Analog Inputs - MER, EEG, EMG.... | | | | |
| 0 | 15 | LFP, EEG (Cont) of Electrode 1..16 | Connector 1, HS 1 on FE Box | Analog Input |
| 16 | 31 | LFP, EEG (Cont) of Electrode 17..32 | Connector 2, HS1 or HS2 on FE Box | Analog Input |
| 32 | 47 | LFP, EEG (Cont) of Electrode 33..48 | Connector 3, HS1 or HS2 on FE Box | Analog Input |
| 48 | 63 | LFP, EEG (Cont) of Electrode 49..64 | Connector 4, HS1 or HS2 on FE Box | Analog Input |
| 64 | 79 | LFP, EEG (Cont) of Electrode 65..80 | Connector 5, HS1 or HS2 on FE Box | Analog Input |
| 80 | 95 | LFP, EEG (Cont) of Electrode 81..96 | Connector 6, HS1 or HS2 on FE Box | Analog Input |
| 96 | 111 | LFP, EEG (Cont) of Electrode 97..112 | Connector 7, HS1 or HS2 on FE Box | Analog Input |
| 112 | 127 | LFP, EEG (Cont) of Electrode 113..128 | Connector 8, HS1 or HS2 on FE Box | Analog Input |
| | | | | |
| 128 | 143 | Segmented of Electrode 1..16 | Connector 1, HS 1 on FE Box | Analog Input |
| 144 | 159 | Segmented of Electrode 17..32 | Connector 2, HS1 or HS2 on FE Box | Analog Input |
| 160 | 175 | Segmented of Electrode 33..48 | Connector 3, HS1 or HS2 on FE Box | Analog Input |
| 176 | 191 | Segmented of Electrode 49..64 | Connector 4, HS1 or HS2 on FE Box | Analog Input |
| 192 | 207 | Segmented of Electrode 65..80 | Connector 5, HS1 or HS2 on FE Box | Analog Input |
| 208 | 223 | Segmented of Electrode 81..96 | Connector 6, HS1 or HS2 on FE Box | Analog Input |
| 224 | 239 | Segmented of Electrode 97..112 | Connector 7, HS1 or HS2 on FE Box | Analog Input |
| 240 | 255 | Segmented of Electrode 113..128 | Connector 8, HS1 or HS2 on FE Box | Analog Input |
| | | | | |
| 256 | 271 | Spike, EMG (Cont) of Electrode 1..16 | Connector 1, HS 1 on FE Box | Analog Input |
| 272 | 287 | Spike, EMG (Cont) of Electrode 17..32 | Connector 2, HS1 or HS2 on FE Box | Analog Input |
| 288 | 303 | Spike, EMG (Cont) of Electrode 33..48 | Connector 3, HS1 or HS2 on FE Box | Analog Input |
| 304 | 319 | Spike, EMG (Cont) of Electrode 49..64 | Connector 4, HS1 or HS2 on FE Box | Analog Input |
| 320 | 335 | Spike, EMG (Cont) of Electrode 65..80 | Connector 5, HS1 or HS2 on FE Box | Analog Input |
| 336 | 351 | Spike, EMG (Cont) of Electrode 81..96 | Connector 6, HS1 or HS2 on FE Box | Analog Input |
| 352 | 367 | Spike, EMG (Cont) of Electrode 97..112 | Connector 7, HS1 or HS2 on FE Box | Analog Input |
| 368 | 383 | Spike, EMG (Cont) of Electrode 113..128 | Connector 8, HS1 or HS2 on FE Box | Analog Input |
| | | | | |
| 384 | 399 | Raw (Cont) of Electrode 1..16 | Connector 1, HS 1 on FE Box | Analog Input |
| 400 | 415 | Raw (Cont) of Electrode 17..32 | Connector 2, HS1 or HS2 on FE Box | Analog Input |
| 416 | 431 | Raw (Cont) of Electrode 33..48 | Connector 3, HS1 or HS2 on FE Box | Analog Input |
| 432 | 447 | Raw (Cont) of Electrode 49..64 | Connector 4, HS1 or HS2 on FE Box | Analog Input |
| 448 | 463 | Raw (Cont) of Electrode 65..80 | Connector 5, HS1 or HS2 on FE Box | Analog Input |
| 464 | 479 | Raw (Cont) of Electrode 81..96 | Connector 6, HS1 or HS2 on FE Box | Analog Input |

| | | | | |
|---|------|-------------------------------------|--|--------------------|
| 480 | 495 | Raw (Cont) of Electrode 97..112 | Connector 7, HS1 or HS2 on FE Box | Analog Input |
| 496 | 511 | Raw (Cont) of Electrode 113..128 | Connector 8, HS1 or HS2 on FE Box | Analog Input |
| | | | | |
| Analog Inputs - General Purpose | | | | |
| 1000 | 1015 | General Purpose Analog Inputs 1..16 | A. In connector on FE Box | Analog Input |
| | | | Analog Inputs Connectors (1-8) on Front Panel of RM Box and duplicated on A. In Connector (Chan 1-8) on Slot 1 on the Rear Panel of the RM Box | |
| 1016 | 1023 | GP Analog Input 17..24 | | Analog Input |
| 1024 | 1031 | GP Analog Input 25..32 | A. In Connector (Chan 9-16) on Slot number 1 on the Rear Panel of the RM Box | Analog Input |
| | | | | |
| User Defined Data Channels | | | | |
| 1032 | 1047 | User Defined Analog Inputs | Internal - Future Use (Movement speed, channel calculation, other) | Analog Input |
| 1050 | 1053 | User Defined Waveform Data | Internal - can be sent to analog outputs or for stimulation | |
| 1060 | 1075 | User Defined Binary Data | Internal - Future Use | |
| | | | | |
| Analog Outputs (Format will be different for scripting purposes) | | | | |
| 1100 | 1107 | GP Analog Output 1-8 | A. Out Connector on Front Panel of FE Box | Analog Output |
| 1108 | 1115 | GP Analog Output 9-16 | Analog Outputs Connectors 1-8 on Front Panel of RM Box | Analog Output |
| | | | | |
| 1120 | 1121 | Audio Output 1, right and left | Audio 1 on Front Panel of FE Box | Analog Output |
| 1122 | 1123 | Audio Output 2, right and left | Audio 2 on Front Panel of FE Box | Analog Output |
| | | | | |
| 1130 | 1131 | Audio Output 1, right and left | Audio 1 on Front Panel of RM Box | Analog Output |
| 1132 | 1133 | Audio Output 2, right and left | Audio 2 on Front Panel of RM Box | Analog Output |
| | | | | |
| | | | | |
| Digital Inputs | | | | |
| 1200 | | Digital Port 1 with stobe and ready | D. In Port Connector on Front Panel of FE Box | Digital Input Port |
| 1201 | | Digital Port 2 | Internal on FE BOX - Not available as a port | |
| 1202 | | Digital Port 3 with stobe and ready | D. In Port Connector on Rear Panel of RM Box - Slot 1 | Digital Input Port |
| 1203 | | Digital Port 4 | Internal on RM BOX - Not available as a port | |
| 1204 | | Digital Port 5 with stobe and ready | D. In Port Connector on Rear Panel of RM Box - Slot 2 | Digital Input Port |

| | | | | |
|------|------|--------------------------------------|--|--------------------|
| 1205 | | Digital Port 6 | Internal on RM BOX - Not available | |
| 1206 | | Digital Port 7 with stobe and ready | D. In Port Connector on Rear Panel of RM Box - Slot 3 | Digital Input Port |
| 1207 | | Digital Port 8 | Internal on RM BOX - Not available | |
| 1208 | | Digital Port 9 with stobe and ready | D. In Port Connector on Rear Panel of RM Box - Slot 4 | Digital Input Port |
| 1209 | | Digital Port10 | Internal on RM BOX - Not available | |
| 1210 | | Digital Port 11 with stobe and ready | D. In Port Connector on Rear Panel of RM Box - Slot 4 | Digital Input Port |
| 1211 | | Digital Port 12 | Internal on RM BOX - Not available | |
| 1212 | | Digital Port 13 with stobe and ready | D. In Port Connector on Rear Panel of RM Box - Slot 6 | Digital Input Port |
| 1213 | | Digital Port 14 | Internal on RM BOX - Not available | |
| 1214 | | Digital Port 15 with stobe and ready | D. In Port Connector on Rear Panel of RM Box - Slot 7 | Digital Input Port |
| 1215 | | Digital Port 16 | Internal on RM BOX - Not available | |
| 1216 | | Digital Port 17 with stobe and ready | D. In Port Connector on Rear Panel of RM Box - Slot 8 | Digital Input Port |
| 1217 | | Digital Port 18 | Internal on RM BOX - Not available | |
| 1218 | | Digital Port 19 with stobe and ready | D. In Port Connector on Rear Panel of RM Box - Slot 9 | Digital Input Port |
| 1219 | | Digital Port 20 | Internal on RM BOX - Not available | |
| 1220 | | Digital Port 21 | Event Data | |
| 1221 | | Digital Port 22 | Stim Marker 1 | |
| 1222 | | Digital Port 23 | Stim Marker 2 | |
| 1223 | | Digital Port 24 | Stim Marker 3 | |
| 1224 | | Digital Port 24 | Stim Marker 4 | |
| | | | | |
| 1300 | 1315 | TTL Input - Port1, bit 1-16 | D.In Port Connector on Front Panel of FE Box | TTL Input |
| 1316 | 1319 | TTL Input - Port 2, bit 1-4 | D.In Connector on Front Panel of FE Box | TTL Input |
| 1320 | 1331 | TTL Input - Port 2, bit 5-16 | Not in use | |
| 1332 | 1347 | TTL Input - Port 3, bit 1-16 | D.In Port Connector on Rear Panel of RM Box - Slot 1 | TTL Input |
| 1348 | 1351 | TTL Input - Port 4, bit 1-4 | Digital Inputs connectors 0-3 on Front Panel of RM Box | TTL Input |
| 1352 | 1363 | TTL Input - Port 4, bit 5-16 | Not in use | |
| 1364 | 1379 | TTL Input - Port 5, bit 1-16 | D.In Port Connector on Rear Panel of RM Box - Slot 2 | TTL Input |
| 1380 | 1395 | TTL Input - Port 6, bit 1-16 | Not in use | |
| 1396 | 1411 | TTL Input - Port 7, bit 1-16 | D.In Port Connector on Rear Panel of RM Box - Slot 3 | TTL Input |
| 1412 | 1427 | TTL Input - Port 8, bit 1-16 | Not in use | |
| 1428 | 1443 | TTL Input - Port 9, bit 1-16 | D.In Port Connector on Rear Panel of RM Box - Slot 4 | TTL Input |
| 1444 | 1459 | TTL Input - Port 10, bit 1-16 | Not in use | |
| 1460 | 1475 | TTL Input - Port 11, bit 1-16 | D.In Port Connector on Rear Panel of RM Box - Slot 5 | TTL Input |
| 1476 | 1491 | TTL Input - Port 12, bit 1-16 | Not in use | |
| 1492 | 1507 | TTL Input - Port 13, bit 1-16 | D.In Port Connector on Rear Panel of RM Box - | TTL Input |

| | | | | |
|------------------------|------|--------------------------------------|---|------------|
| | | | Slot 6 | |
| 1508 | 1523 | TTL Input - Port 14, bit 1-16 | Not in use | |
| 1524 | 1539 | TTL Input - Port 15, bit 1-16 | D.In Port Connector on Rear Panel of RM Box - Slot 7 | TTL Input |
| 1540 | 1555 | TTL Input - Port 16, bit 1-16 | Not in use | |
| 1556 | 1571 | TTL Input - Port 17, bit 1-16 | D.In Port Connector on Rear Panel of RM Box - Slot 8 | TTL Input |
| 1572 | 1587 | TTL Input - Port 18, bit 1-16 | Not in use | |
| 1588 | 1603 | TTL Input - Port 19, bit 1-16 | D.In Port Connector on Rear Panel of RM Box - Slot 9 | TTL Input |
| 1604 | 1619 | TTL Input - Port 20, bit 1-16 | Not in use | |
| | | | | |
| Digital Outputs | | | | |
| 1700 | | GP Digital Output Port 1 (Bits 1-8) | D. Out Connector on Front Panel of FE Box | TTL Output |
| 1701 | | GP Digital Output Port 2 (Bits 1-8) | Digital Outputs Connectors 1-8 on Front Panel of RM Box | TTL Output |
| 1702 | | GP Digital Output Port 3 (Bits 1-8) | Internal on ADIO Board installed in Slot 2 of RM Box | TTL Output |
| 1703 | | GP Digital Output Port 4 (Bits 1-8) | Internal on ADIO Board installed in Slot 3 of RM Box | TTL Output |
| 1704 | | GP Digital Output Port 5 (Bits 1-8) | Internal on ADIO Board installed in Slot 4 of RM Box | TTL Output |
| 1705 | | GP Digital Output Port 6 (Bits 1-8) | Internal on ADIO Board installed in Slot 5 of RM Box | TTL Output |
| 1706 | | GP Digital Output Port 7 (Bits 1-8) | Internal on ADIO Board installed in Slot 6 of RM Box | TTL Output |
| 1707 | | GP Digital Output Port 8 (Bits 1-8) | Internal on ADIO Board installed in Slot 7 of RM Box | TTL Output |
| 1708 | | GP Digital Output Port 9 (Bits 1-8) | Internal on ADIO Board installed in Slot 8 of RM Box | TTL Output |
| 1709 | | GP Digital Output Port 10 (Bits 1-8) | Internal on ADIO Board installed in Slot 9 of RM Box | TTL Output |
| | | | | |
| 1800 | | Flushing channel | Internal Use | |
| | | | | |
| 1999 | | Steam Data | Internal Use | |
| | | | | |